

RESUSCITAIRE CHECKLIST WRH A&E

Top Surface	
Heater Light Clock Clean, flat working surface Power Supply O2 and Air supply with pressure relief valve (replace if < 500 psi) Suction apparatus with yankauer	DAILY CHECK TO ENSURE WORKING
Maternity delivery pack	2
Warm towels	2
Hats – large & small	1 of each
Stethoscope	1
Breathing Circuit with T-Piece	1
Masks – size 00, 0.5, 0, 1	1 of each
Self-inflating BVM bag (~500ml) with blow off valve	1 (hanging on side of resus)
Oxygen Saturation probe	1
Cord clamp	2
Top Drawer	
Guedel Airway – size 1, 0, 00 & 000	1 of each
Laryngoscope Handle	1
Laryngoscope Blades – size 00, 0, 1	1 of each
Tendertip suction catheters- 5Fr, 6Fr, 8Fr, 10Fr	2 of each
Endotracheal Intubation Stylet	2
Endotracheal Tubes – size 2.0, 2.5, 3.0, 3.5	2 of each
End-tidal CO2 Sensor	1
Medi-dermaS barrier applicator	1 box
Neofit ETT fixation	2
i-Gel size 1	1
AirQ3 LMA size 0	1
Optilube sachets	2
Surfactant Administration Kit	1
Bottom Drawer	
Box non-sterile gloves	1
Spare towels	2
Spare hats – large/ small	1 of each
Spare neonatal face masks- 00, 0, 0.5, 1	1 of each
Spare yankauer	1
Spare breathing circuit with T-piece	1
Neohelp bags- <1kg and <2kgs	1 of each
Digital Thermometer and probe covers	1
Baby labels	2
Other	
Checklist folder with signing sheets	
Neonatal algorithms and scribe sheets	

Blue Emergency Drug Box Kept In Resus Cupboard



RESUSCITAIRE CHECKLIST WRH A&E

