

RESUSCITAIRE CHECKLIST

| Top Surface | |
|---|----------------------------------|
| Heater Light Clock Clean, flat working surface Power Supply O2 and Air supply with pressure relief valve (replace if < 500 psi) Suction apparatus with yankauer | DAILY CHECK TO ENSURE WORKING |
| Warm towels | 2 |
| Hats – large & small | 1 of each |
| Stethoscope | 1 |
| Breathing Circuit with T-Piece | 1 |
| Masks – size 00, 0.5, 0, 1 | 1 of each |
| Self-inflating BVM bag (~500ml) with blow off valve | 1 (hanging on side of resus) |
| Oxygen Saturation probe | 1 |
| Top Drawer | |
| Neonatal face masks- 00, 0, 0.5, 1 | 1 of each |
| Laryngoscope Handle | 2 |
| Laryngoscope Blades – size 00, 0, 1 | 1 of each |
| Yankauer | 2 |
| Tendertip suction catheters- 5Fr, 6Fr, 8Fr, 10Fr | 2 of each |
| Endotracheal Intubation Stylet | 2 |
| Endotracheal Tubes – size 2.0, 2.5, 3.0, 3.5 | 2 of each |
| End-tidal CO2 Sensor | 2 |
| Medi-dermaS barrier applicator | 1 box |
| Neofit ETT fixation | 2 |
| i-Gel size 1 | 1 |
| AirQ3 LMA size 0 | 1 |
| Optilube sachets | 2 |
| Surfactant Administration Kit | 2 |
| Bottom Drawer | |
| Box non-sterile gloves | 1 |
| Spare towels | 2 |
| Spare hats – large/ small | 1 of each |
| Spare breathing circuit with T-piece | 1 |
| Neohelp bags- <1kg and <2kgs | 1 of each |
| Digital Thermometer and probe covers | - |
| Other | |
| Checklist folder with neonatal algorithms | |
| Scribe clipboard and sheets | |
| Baby Labels | |

Blue Emergency Drug Box Kept

